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**INSTRUCTIONS FOR USE**  
**Bimanual I/A System Disposable – EN**

**CONTENT**

Each set contains two sterile instruments (handles). The handles are designed with standard Luer-cones. According to the type of Luer-cone (female/male) and the colour code, the handles can be used for irrigation and aspiration, respectively (see table):

Set content	20G / 21G	22G / 23G / 25G
Irrigation Handle (female Luer)	Blue	Violet
Aspiration Handle (male Luer)	Red	Transparent

**INTENDED USE**

The Bimanual I/A Disposable Instruments are intended to be used to deliver solutions to the eye and to extract/aspirate solutions and residuals of tissue from the eye, respectively, during cataract operations in the anterior chamber. The instruments are connected to the aspiration / irrigation system by standard Luer-cones. The Bimanual I/A Disposable Instruments are intended to be used only by appropriately specialist surgeon.

**IMPORTANT USER INFORMATION**

Insufficient flow conditions can cause an imbalance of the fluids. It is strongly recommended to observe the pressure conditions and, if necessary, to adjust them. It is very important to prevent any collapse of the anterior chamber. Bürki inno med AG disclaims any liability for any inappropriate and incorrect handling of the sterile bimanual I/A instruments.

**CHECK BEFORE USE**

Following inspection are necessary before use:

- Visual inspection of obstruction, deformation and no cracks
- Filling of irrigation and aspiration system with NaCl 0.9% solution to ensure the functionality

**WARNINGS/ PRECAUTIONS** (see also packaging symbols)

- The Bimanual I/A Disposable Instruments are designed for single use; any **re-sterilisation or re-use is strongly forbidden and can lead to dangerous post-operative infections of the eye.**
- Don't use the instruments after expiry date or if the packaging is damaged.
- Prior to use, inspect the tips of the Bimanual I/A System. Don't use the product, if the tips are damaged.
- The instruments are only to be used by qualified professionals.
- The manufacturer shall not be liable for any injury or damage suffered by a patient due to the use of the product.
- The Bimanual I/A System are packed in a blister and a Tyvek bag. The Tyvek bag represents the sterile barrier. The blister is an additional protection which is intended to facilitate the handling in the operating room.
- **BSD 600**; as elliptical system it has a higher impermeability on the edge of the cannula. This leads to substantial less efflux. In consequence, specially the **irrigation** has the potential of a substantially increase of the inner eye pressure due to the specific TWIN PORT design. To prevent any damage, e.g. any form of collapse of the eye, a reduced height of the irrigation liquid bottle is strongly recommended. Keep the pressure conditions in permanent observation and, if necessary, adjust (decrease).

**TECHNICAL SPECIFICATION**

Parameter	Specification
Life span	< 60min
Connectors	Standard Luer-cones
Shelf life	3 years in original packaging
Sterilisation type	Sterilization by irradiation
State of delivery	Sterile
Quality Management System of Bürki inno med AG	EN ISO 13485

**DISPOSAL**

After usage, the Bimanual I/A Disposable Instruments must be handled as a contaminated item.

The devices must be disposed according to local regulations. For information, contact your local environmental or public office and appropriate waste disposal companies. In this regard, waste is to be recycled or disposed of:

- without danger to human health
- without using procedures or methods harmful to the environment, particularly to the water, air, soil, flora and fauna
- without creating noise or smell

**DEVICE LABEL AND SYMBOLS**

	Single use only	<b>REF</b>	Article number		Do not resterilize
	Do not use after expiry date	<b>STERILE R</b>	Sterilisation by irradiation		Flow direction
<b>LOT</b>	Batch / Lot Number		Consult instructions for use		Manufacturer
<b>EC REP</b>	Authorised representative Bürki inno med GmbH Im Schaffner 47/1, 69123 Heidelberg-Pfaffengrund, Germany		Do not use if the packaging has been damaged		ID number of notified body
	Single sterile barrier system	<b>Manufacturer:</b>	Bürki inno med AG · Industriestrasse 67 · CH-9443 Widnau · Switzerland Tel. +41 (0)71 744 88 66 · www.buerki.eu		

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Bürki inno med AG rejette toute responsabilité pour les applications inappropriées des instruments bimanuels.