**GEBRAUCHSANWEISUNG**  
**Coaxial I/A System Disposable - DE**

**INHALT**

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**INSTRUCTIONS FOR USE**  
**Coaxial I/A System Disposable - EN**

**CONTENT**

Each bag contains a sterile, coaxial instrument for irrigation and aspiration.

**INTENDED USE**

The Coaxial (mono-manual) I/A Disposable Instrument is intended to be used to deliver solutions to the eye and to extract/aspirate solutions and residuals of tissue from the eye, respectively, during cataract operations in the anterior chamber. The instrument is connected to the aspiration / irrigation system by two standard Luer-cones. The Coaxial (mono-manual) I/A Disposable Instrument is intended to be used only by appropriately specialised surgeons.

**IMPORTANT USER INFORMATION**

Insufficient flow conditions can cause an imbalance of the fluids. It is strongly recommended to observe the pressure conditions and, if necessary, to adjust them. It is very important to prevent any collapse of the anterior chamber. Bürki inno med AG disclaims any liability for any inappropriate and incorrect handling of the sterile coaxial I/A instrument.

**CHECK BEFORE USE**

Following inspections are necessary before use:

- Visual inspection of obstruction, deformation and small cracks
- Fill the irrigation and aspiration system with NaCl 0.9% solution to ensure the functionality

**WARNINGS/ PRECAUTIONS** (see also the packaging symbols)

- The Coaxial (mono-manual) I/A Disposable Instrument is designed for single use; any **re-sterilisation or re-use is strongly forbidden and can lead to dangerous post-operative infections.**
- Don't use the instrument after expiry date or if the packaging is damaged.
- Prior the use, inspect the tip of the Coaxial (mono-manual) I/A System. Don't use the product, if the tip is damaged.
- The instrument is only to be used by qualified professionals.
- The manufacturer is not liable for any injury or damage caused to a patient the use of the product.
- The coaxial I/A System is packed in a Blister and a Tyvek-Bag. The Tyvek bag represents the sterile barrier. The blister is an additional protection which is intended to facilitate the handling in the operating room.

**TECHNICAL SPECIFICATIONS**

Parameter	Specification
Duration of use	< 60min
Connectors	Standard Luer-cones
Shelf life	3 years in original packaging
Sterilisation type	Sterilization by irradiation
State of delivery	Sterile
Quality Management System of Bürki inno med AG	EN ISO 13485

**DISPOSAL**

After usage, the Coaxial I/A Disposable Instrument must be handled as a contaminated item.

The device must be disposed according to local regulations. For information, contact your local environmental or public office and appropriate waste disposal companies. In this regard, waste is to be recycled or disposed of:

- without danger to human health
- without using procedures or methods harmful to the environment, particularly to the water, air, soil, flora and fauna
- without creating noise or smell

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Bürki inno med AG rejette toute responsabilité pour les applications inappropriées de l'instrument COAXIAL stérile.