

Bimanual I/A System reusable: Angled Aspiration

- Lightweight titanium handpieces
- Compatible with Phaco or I/A System via standard Luer connectors
- Angled aspiration handpieces left and right available
- 25G
- Swiss Design / Swiss Made



Tip variants



Angled Aspiration Handpiece left (L)



Angled Aspiration Handpiece right (R)



Angled aspiration handpieces

Art. No.	Description
BA 260/A	Aspiration Handpiece right (R) 25G (Ø0.50mm), asp. port 0.60x0.25mm
BA 261/A	Aspiration Handpiece left (L) 25G (Ø0.50mm), asp. port 0.60x0.25mm

Please note: The expected lifetime of our reusable medical devices is validated by a defined number of processing cycles. Please refer to the information provided in the instructions for use.