

## Bimanual I/A System single-use: Angled Aspiration

- Lightweight handpieces
- Compatible with Phaco or I/A System via standard Luer connectors
- Angled aspiration handpieces left and right available
- 25G
- Swiss Design / Swiss Made



### Tip variants



Angled Aspiration Handpiece left (L)



Angled Aspiration Handpiece right (R)



### Angled aspiration handpieces

Art. No.	Description
BSD 530	Bimanual I/A System single-use 25G (Ø0.50mm), 1 aspiration handpiece angled right: asp. port 0.60x0.25mm, 1 aspiration handpiece angled left: asp. port 0.60x0.25mm