

# Monomanual I/A System single-use (coaxial)

- Lightweight monomanual irrigation / aspiration handpiece
- Compatible with Phaco or I/A System via standard Luer connectors
- Smooth abrasive (rough) tip available for capsule polishing
- Angled tip: 45° or 65° available



## Tip variants



45° · angled



45° · Soft Tip



65° · angled

## Monomanual I/A System single-use, angled 45°

| Art. No.        | Description                                                                                                                                                                              |
|-----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>CX2045-P</b> | Monomanual I/A System single-use 20G (Ø0.90mm), polished tip, angled 45°, asp. port Ø0.35mm, incl. silicone sleeve 16G (Ø1.70mm) for irrigation with dual-side ports Ø1.10mm             |
| <b>CX2045-R</b> | Monomanual I/A System single-use 20G (Ø0.90mm), rough tip, angled 45°, asp. port Ø0.35mm, incl. silicone sleeve 16G (Ø1.70mm) for irrigation with dual-side ports Ø1.10mm                |
| <b>CX2145-P</b> | Monomanual I/A System single-use 21G (Ø0.80mm), polished tip, angled 45°, asp. port Ø0.35mm, incl. silicone sleeve 17G (Ø1.50mm) for irrigation with dual-side ports 1.65x0.77mm         |
| <b>CX2145-R</b> | Monomanual I/A System single-use 21G (Ø0.80mm), rough tip, angled 45°, asp. port Ø0.35mm, incl. silicone sleeve 17G (Ø1.50mm) for irrigation with dual-side ports 1.65x0.77mm            |
| <b>CX2145-S</b> | Monomanual I/A System single-use 21G, Soft Tip System, angled 45°, incl. silicone sleeve 17G (Ø1.50mm) for irrigation with dual-side ports 1.60x0.83mm and aspiration, asp. port Ø0.38mm |
| <b>CX2245-P</b> | Monomanual I/A System single-use 22G (Ø0.70mm), polished tip, angled 45°, asp. port Ø0.35mm, incl. silicone sleeve 18G (Ø1.40mm) for irrigation with dual-side ports 1.50x0.67mm         |
| <b>CX2245-R</b> | Monomanual I/A System single-use 22G (Ø0.70mm), rough tip, angled 45°, asp. port Ø0.35mm, incl. silicone sleeve 18G (Ø1.40mm) for irrigation with dual-side ports 1.50x0.67mm            |

# Monomanual I/A System single-use (coaxial)

**Monomanual I/A System single-use, angled 65°**

| Art. No.        | Description                                                                                                                                                                      |
|-----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>CX2065-P</b> | Monomanual I/A System single-use 20G (Ø0.90mm), polished tip, angled 65°, asp. port Ø0.35mm, incl. silicone sleeve 16G (Ø1.70mm) for irrigation with dual-side ports Ø1.10mm     |
| <b>CX2065-R</b> | Monomanual I/A System single-use 20G (Ø0.90mm), rough tip, angled 65°, asp. port Ø0.35mm, incl. silicone sleeve 16G (Ø1.70mm) for irrigation with dual-side ports Ø1.10mm        |
| <b>CX2165-P</b> | Monomanual I/A System single-use 21G (Ø0.80mm), polished tip, angled 65°, asp. port Ø0.35mm, incl. silicone sleeve 17G (Ø1.50mm) for irrigation with dual-side ports 1.65x0.77mm |
| <b>CX2165-R</b> | Monomanual I/A System single-use 21G (Ø0.80mm), rough tip, angled 65°, asp. port Ø0.35mm, incl. silicone sleeve 17G (Ø1.50mm) for irrigation with dual-side ports 1.65x0.77mm    |

