

## ANGLED Aspiration Instruments reusable and single use 25G



### Product Features reusable

- Titanium lightweight, colour-coded Handpieces
- Standard male/female Luer connector
- Compatible with Phaco or I/A System
- Smooth bullet shaped tip
- 25 G

### Product Features single use

- colour-coded Handpieces
- Standard male/female Luer connector
- Compatible with Phaco or I/A System
- Smooth bullet shaped tip
- Pack of left and right disposable aspiration handle
- 25 G

### Aspiration

#### BA 260/A

1 Aspiration Handle right (R),  
cannula 25G (0.5mm), Asp. Port 0.25x0.6mm


#### BA 261/A

1 Aspiration Handle left (L),  
cannula 25G (0.5mm), Asp. Port 0.25x0.6mm

### Aspiration

#### BSD 530

Aspiration Handle right (R), cannula 25G (0.5mm),  
Asp. Port 0.25x0.6mm  
Aspiration Handle left (L), cannula 25G (0.5mm),  
Asp. Port 0.25x0.6mm



The expected lifetime of our reusable medical devices is validated by a defined number of processing cycles. Please refer to the information provided in the instructions for use.

Version 11/2024

 Swiss Design / Swiss Made

About Bürki  
Bürki plans, develops, produces microsurgical  
instruments in SWISS QUALITY and according to  
EN ISO 13485

Bürki inno med AG  
Industriestrasse 67  
9443 Widnau  
Switzerland

Tel +41 (0)71 744 88 66  
info@buerki.eu  
www.buerki.eu