

## Bimanual System Disposable Irrigation Aspiration 22G, 23G



### 22G



#### **BSD 540**

Bimanual System Disposable **22G** consisting of:

1 Irrigation Handpiece, with curved cannula and dual side ports  
ø 0.45mm

1 Aspiration Handpiece with curved cannula and asp. port ø 0.35 mm,  
smooth abrasive (rough) tip



#### **BSD 541**

Bimanual System Disposable **22G** consisting of:

1 Irrigation Handpiece, with curved cannula and dual side ports  
ø 0.45 mm

1 Aspiration Handpiece with curved cannula and asp. port ø 0.35 mm,  
polished tip



#### **BSD 542**

Bimanual System Disposable **22G** consisting of:

1 Irrigation Handpiece, with curved cannula and open end flushing 45°  
polished

1 Aspiration Handpiece with curved cannula and asp. port ø 0.35 mm,  
smooth abrasive (rough) tip

### 23G



#### **BSD 522**

Bimanual System Disposable **23G** consisting of:

1 Irrigation Handpiece with curved cannula and open end flushing 45°,  
polished

1 Aspiration Handpiece with curved cannula and asp. port ø 0.30 mm,  
smooth abrasive (rough) tip

### ***BSD 520 and BSD 521 now available***



#### **BSD 520**

Bimanual System Disposable **23G** consisting of:

1 Irrigation Handpiece, with curved cannula, dual side ports ø 0.35 mm

1 Aspiration Handpiece, with curved cannula, Asp. Port ø 0.30 mm  
smooth abrasive (rough) tip



#### **BSD 521**

Bimanual System Disposable **23G** consisting of:

1 Irrigation Handpiece with curved cannula, dual side ports ø 0.35 mm

1 Aspiration Handpiece with curved cannula, Asp. Port ø 0.30 polished  
tip

**Detailed product Features:**

	<b>BSD 540</b>	<b>BSD 541</b>	<b>BSD 542</b>	<b>BSD 520</b>	<b>BSD 521</b>	<b>BSD 522</b>
Single use / no risk of cross contamination	●	●	●	●	●	●
Lightweight, colour-coded handpiece	●	●	●	●	●	●
Standard male/female Luer connector	●	●	●	●	●	●
Compatible with your Phaco or I/A System	●	●	●	●	●	●
Smooth bullet-shaped tip		●			●	
Irrigation front flushing 45°			●			●
Irrigation dual side ports 0.35 mm				●	●	
Irrigation dual side ports 0.45 mm	●	●				
Aspiration port 0.30 mm, rough tip				●		●
Aspiration port 0.30 mm, smooth tip					●	
Aspiration port 0.35 mm, rough tip	●		●			
Aspiration port 0.35 mm, smooth tip		●				

Version 08/2024

About Bürki  
 Bürki plans, develops, produces and certifies  
 microsurgical instruments in SWISS QUALITY  
 DIN EN ISO 13485

Bürki inno med AG  
 Industriestrasse 67  
 9443 Widnau  
 Switzerland

Tel +41 (0)71 744 88 66  
 Fax +41 (0)71 744 88 67  
 info@buerki.eu  
 www.buerki.eu