

## Reusable Bimanual Irrigation Aspiration System 22G, 23G



### Product Features:

- Lightweight, colour coded irrigation and aspiration Handpieces
- Compatible with your existing Phaco or I/A System
- Smooth abrasive (rough) tip available for capsule polishing
- Swiss Design / Swiss Made

### Aspiration

#### 22G




**BA 118/T**  
Aspiration Handpiece, polished tip,  
(0.70mm), Asp. port 0.35 mm



**BA 119/T Special model**  
Aspiration Handpiece, smooth  
abrasive (rough) tip, (0.70mm)  
Asp. port 0.35 mm

#### 23G



**BA 113/T**  
Aspiration Handpiece, polished tip,  
(0.65mm) Asp. port 0.30 mm



**BA 114/T**  
Aspiration Handpiece, smooth  
abrasive (rough) tip, (0.65mm)  
Asp. port 0.30 mm

### Irrigation



**BI 128/T Special model**  
Irrigation Handpiece, polished tip,  
dual side ports 0.45 mm



**BI 129/T Special model**  
Irrigation Handpiece, polished tip,  
open end Tip 45°



**BI 123/T**  
Irrigation Handpiece, round cannula,  
polished tip, dual 0.35 mm side ports



**BI 124/T**  
Irrigation Handpiece, polished tip  
open end Tip 45° (0.65mm)

The expected lifetime of our reusable medical devices is validated by a defined number of processing cycles. Please refer to the information provided in the instructions for use.

Version 11/2024

About Bürki  
Bürki plans, develops, produces and certifies  
microsurgical instruments in SWISS QUALITY  
EN ISO 13485

Bürki inno med AG  
Industriestrasse 67  
9443 Widnau  
Switzerland

Tel +41 (0)71 744 88 66  
Fax +41 (0)71 744 88 67  
[info@buerki.eu](mailto:info@buerki.eu)  
[www.buerki.eu](http://www.buerki.eu)