

Bimanual Reusable Family Irrigation 20G, 21G, 22G, 23G

Product Features:

- Lightweight, colour coded irrigation Handpieces
- Compatible with your existing Phaco or I/A System
- Smooth abrasive (rough) tip available for capsule polishing
- Swiss Design / Swiss Made

Irrigation



20G



BI 126/T
Irrigation Handle open end tip
45° (0.90mm), High flow

21G



BI 120/T
Irrigation Handle polished tip
(0.80mm), dual side ports
0.50mm



BI 125/T
Irrigation Handle smooth
abrasive (rough) tip,
(0.80mm), dual side ports
0.50mm



BI 140/T
Irrigation Handle polished tip
open end 45°

22G



BI 128/T special model
Irrigation Handle polished tip,
dual side ports 0.45mm



BI 129/T special model
Irrigation Handle polished tip,
open end 45°

23G



BI 123/T
Irrigation Handle polished tip
(0.65mm), dual side ports
0.35mm



BI 124/T
Irrigation Handle polished tip
(0.65mm), open end 45°



BI 127/T special model (MLI)
Irrigation Handle polished tip
(0.65mm)

The expected lifetime of our reusable medical devices is validated by a defined number of processing cycles. Please refer to the information provided in the instructions for use.

Version 11/2024