

Bimanual System Disposable Irrigation Aspiration 22G, 23G



22G



BSD 540

Bimanual System Disposable **22G** consisting of:

1 Irrigation Handpiece, with curved cannula and dual side ports
ø 0.45mm

1 Aspiration Handpiece with curved cannula and asp. port ø 0.35 mm,
smooth abrasive (rough) tip



BSD 541

Bimanual System Disposable **22G** consisting of:

1 Irrigation Handpiece, with curved cannula and dual side ports
ø 0.45 mm

1 Aspiration Handpiece with curved cannula and asp. port ø 0.35 mm,
polished tip



BSD 542

Bimanual System Disposable **22G** consisting of:

1 Irrigation Handpiece, with curved cannula and open end flushing 45°
polished

1 Aspiration Handpiece with curved cannula and aps. port ø 0.35 mm,
smooth abrasive (rough) tip

23G



BSD 522

Bimanual System Disposable **23G** consisting of:

1 Irrigation Handpiece with curved cannula and open end flushing 45°,
polished

1 Aspiration Handpiece with curved cannula and aps. port ø 0.30 mm,
smooth abrasive (rough) tip

Detailed product Features:

	BSD 540	BSD 541	BSD 542	BSD 522
Single use / no risk of cross contamination	●	●	●	●
Lightweight, colour-coded handpiece	●	●	●	●
Standard male/female Luer connector	●	●	●	●
Compatible with your Phaco or I/A System	●	●	●	●
Smooth bullet-shaped tip		●		
Irrigation front flushing 45°			●	●
Irrigation dual side ports 0.35 mm				
Irrigation dual side ports 0.45 mm	●	●		
Aspiration port 0.30 mm, rough tip				●
Aspiration port 0.30 mm, smooth tip				
Aspiration port 0.35 mm, rough tip	●		●	
Aspiration port 0.35 mm, smooth tip		●		

Version 12/2024

About Bürki
 Bürki plans, develops, produces and certifies
 microsurgical instruments in SWISS QUALITY
 DIN EN ISO 13485

Bürki inno med AG
 Industriestrasse 67
 9443 Widnau
 Switzerland

Tel +41 (0)71 744 88 66
 Fax +41 (0)71 744 88 67
 info@buerki.eu
 www.buerki.eu